

LED Light Source 24Watt



Accepts standard Schott—Fostec light guides.

Ferrule size is 0.718" (5/8")

No fan.

ON/OFF switch with dimmer control.

Smaller footprint; (3.75" W x 6" D x 5.75" H)

Compatible with most fiber optic cables or guides

Life: 20,000 hrs

Color temperature: 6000~7000k

LED power: 24W

Filters available: Yellow, Red, Blue

Power input 96~265 VAC

Light guides include:

- Single flexible,
- Single gooseneck,
- Twin gooseneck
- Continuous Ring
- Continuous Ring (40mm id)

imageoptics.co.uk

www.close-ups.co.uk

Tel 00 44(0)118 9793727