

Track Stands



Tall Track Stand

76mm Scope Holder for microscopes or lenses

Vertical Column: 500 mm

Travel Distance: 375 mm

Fine/Course Focus mechanism: Fine 0.56 mm per rotation,
Course 23 mm per rotation

LED Base: 320 x 316 x 21 mm

95 mm dia Glass Plate & 95 mm black/white contrast plate

90—240v 50/60 Hz supply



Track Stand 45 mm holder

Base: 310 w x 305 d x 16mm

Column height: 235 mm

Travel: 140 mm

Scope holder: 45mm diameter

Centre holder to centre of column 140 mm



With 6000 Macrozoom lens (adapter required) and
mechanical XY stage.

Shown with co-axial light option for the inspection
of semiconductors and cross-section.

Also compatible with 12x Zoom (Macro or Micro),
LED Ring lights, co-axial lights and measuring stage
with 25 mm digital micrometers.