

VM12x Video Microscope

A Quality Zoom system for the most demanding applications

The **VM12x** offers the highest combination of zoom range and resolution in an optical inspection system. It is Ideal for medium to high magnification and gives a large zoom range, increased resolution and unmatched edge flatness and clarity.

High Quality Optics for Research applications



A modular zoom system for:

High Speed Cameras


High Magnification

Variable working distances

White light, IR or Polarising

Inspection or Measurement

MEMS, Microfluidics



- Image capture
- Video Capture
- High zoom ration
- High Resolution HDMI camera
- Measurements on screen
- Exported Excel file
- No computer required for operation

Stands:

Boom stands, bench stands and horizontal viewing stand with fine focus XYZ

Options:

LED Ring Lights + LED backlights and side lights

Detents in the zoom to aid accurate measurements

Adapters: 0.5x, 0.67x, 1x, 2x , 3.5x. **(C + F-mount options)**

Lens attachments: 0.25x, 0.5x, none, 0.75x, 1.5x and 2x.