


Video Measurement system



Larger 1/2" Sony sensor - Greater sensitivity and larger field of view

Auto calibration

Measure without a computer

Full HD 1980 x 1080 Colour HDMI camera

1/2" Sony sensor for image quality and high sensitivity

Export excel spreadsheet and **image with measurements**

Menus: camera controls, line width/colour

Digital zoom x4, x9, x16 to aid cursor placement

Wireless mouse, psu, cable and calibration slide included.

Suitable for all C-mount lenses (including telecentric)

Suitable for microscopes (stereo and compound)

See... imageoptics.co.uk for sample captured images

Take multiple measurements in different layers.

Move, amend, delete individual measurements

Annotate images

Live HD images to monitor with no delay.

Multiple calibrations

Measurements:

Mark and annotation

Horizontal line

Vertical Lines

Point to point

Rectangular

Circle


Fan shaped (three point)

Circle centre to centre

Concentric circles

Angle measurement

Polygon



Our new website

imageoptics.co.uk

enquiries@imageoptics.co.uk

74 Reading Rd, Wokingham, Berkshire, RG41 1EL

Telephone 0118 979 3727

Details may change