

# LED illumination



## LED RING Light for Stereo Microscopes and Macro Zoom Lenses


60 White LED's, 6000K

60 mm internal diameter with three fixing screws.

ABS enclosure, dimmer is also the off switch

96 to 220 V **£130**

Adapters for 10:1 zoom, Meiji EMT1 or 2, Navitar 12x **£28**



## 4 Zone control Ring Light

60 White LED's, 6400K

60 mm internal diameter with three fixing screws

Working distance 50 mm to 150 mm

On/Off, dimmer and quadrant control on ring light

Includes 12v dc psu **£145**



## LED Ring light for smaller lenses

A perfect fit for smaller lenses such as the 3:1 zoom or fixed focal lenses (eg 12mm 2/3" )

40 White LED's

Working Distance = 40 mm to 200 mm

95 to 240v , 5 watt, 6000 K, ABS enclosure **£ 130**



## Twin LED illuminator

A twin gooseneck LED illuminator. Each arm has a 3 watt white LED and focus lens. The arms are 590 mm long and there is on/off and dimmer to give fine control. Supplied with UK mains adapter. Colour temp. 3200 degree kelvin. **£145**



## 8 point LED ring Light

8 Watt, 6000 to 7000 K

UV Free, ESD Safe

Inner diameter 60 mm

96 to 265 V , 30 vdc

50,000 hour life **£210**

Optional polarising set  
for 8 point ring light

**£85**



Our new website

[imageoptics.co.uk](http://imageoptics.co.uk)

[enquiries@imageoptics.co.uk](mailto:enquiries@imageoptics.co.uk)

74 Reading Rd, Wokingham, Berkshire, RG41 1EL

Telephone 0118 979 3727

Copyright © Close-Ups 2017 prices exclude VAT. Prices may vary