

Research lenses & cameras

Best prices for Navitar, Mitutoyo, Lumenera, IDS, Computar, Kowa, Photron



Navitar MTL

This **NEW** modular tube lens system offers the ideal optical solution for OEM, research applications, metrology, flat panel inspection and cell imaging.

The MTL can be used with **C-mount, F-mount, SLR or 42 mm** mount cameras.

The lens system can be mounted in bench stands, in boom stands, on optical tables etc.


The MTL is ideal for Fluorescent imaging, High Speed video or co-axial illumination applications. Also a range of LED illumination.



Navitar 12x and 6000 (zoom) and Precise Eye (non-zoom)

The 12x Zoom with its every growing range of adapters and accessories has proven to be a successful optical tool in many manufacturing and research centres.

- **High Speed imaging in both Macro and Micro applications**
- **With C or F mount cameras**
- **Macro and Micro applications with cameras up to 1" format**
- **Extended working distance options (beyond a metre)**
- **Near UV and Near IR options - lenses and Objectives**



Photron High Speed Cameras

Fastcam mini UX 50—2,000 fps @ 1028 x 1024 to 160,000 fps @ 1280 x 8

Fastcam mini UX 100—4,000 fps @ 1280 x 1024 to 800,000 fps @ 640 x 8



Lumenera
Microscopy
Cameras

| Model | | (MegaPixel) | Sensor | Sensor size | Interface |
|-----------------|--------|-------------|----------------|-------------|-------------|
| Infinity 3-3UR | C or M | 2.8 | icx 674 Had II | 2/3" | USB3 |
| Infinity 3-3URF | C or M | 3 | icx 675 Had II | 2/3" | USB3 |
| Infinity 3-6UR | C or M | 6 | icx 694 Had II | 1" | USB3 |
| Infinity 3S-1UR | C or M | 1.4 | icx 825 Had II | 2/3" | USB3 |



IDS
Microscopy /
Machine Vision
Cameras

| | | | | | | |
|-----------|--------|-------------------|------------|---------|--------|---------------|
| UI-3270LE | C or M | 3.2 | imx 265 | 3.45 um | 1/1.8" | USB3.1 |
| UI-3260CP | C or M | 2.4 | imx 249 | 5.86 um | 1/1.2" | USB3 |
| UI-3280CP | C or M | 5 | imx 264 | 3.45 um | 2/3" | USB3 |
| UI-3130CP | C or M | 800x600 (575 fps) | Python 500 | 4.8 um | 1/3.6" | USB3 |



Our new website imageoptics.co.uk

Wokingham, Berkshire, RG4 1EL

enquiries@imageoptics.co.uk

Tel +44 (0) 118 979 3727

Established 1994 Copyright © close-ups 2019