

CU4K8M HDMI Colour Camera



NEW High Resolution camera
for lenses and microscopes

C-mount compatible for popular lenses and microscopes.

With a high resolution and high sensitivity Starvis Sony sensor the **NEW** CU4K8M camera offers live imaging and image and video capture directly to USB drive.

A user friendly on-screen menu allows browsing captured images, image comparison, line generation and image measurement. Image measurements can be exported in Excel format.

A wired or wireless mouse is connected to the camera to access the user interface and a second USB port accepts a USB flash drive.

HDMI specification	
SENSOR	IMX334
Sensor size	1/1.8
Pixel size	2.0*2.0µm
Output resolution	3840*2160
Frame rate	HDMI+USB@ 30FPS
Image storage	USB stick or computer
Video	1920*1080&3840*2160@30FPS
Exposure time	0-999MS
Exposure	Auto and manual, providing target brightness settings
White balance	Auto and manual
Parameter	Contrast, saturation, gamma, sharpness, noise reduction etc.
Working temperature	0-60°C
Storage temperature	-20-85°C
Size	mm
N.W.	0.5kg
Power	DC 12V

imageoptics.co.uk Telephone 01189793727

Wokingham, Berkshire RG41 1EL